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.....  
Typ stetoskopu

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Sprzedawca

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### **SPL™ STETHOSCOPES** Instructions for use

The product is a medical diagnostic tool designed for auscultation of the patient. Read these instruction carefully before using the stethoscope and keep them in a safe place for future reference.

#### Setting-up

Plastic parts don't consist latex. The stethoscope is supplied with soft eartips fitted. These can be unscrewed and replaced. Stethoscope binaural should be installed so eartip body work around axis aural. Eartips bodies regulation pressure gets by binaural arms extension or squeezing them within. The chestpiece lets you choose between a membrane and an open cone. To select, hold the tube in one hand and rotate the chestpiece with the other hand 180° until it clicks into place. Position the chestpiece without applying too much pressure.

#### Cleaning and maintenance

Eartips bodies and acoustic tubing should be often washed in hot water with soap addition towards removing fat from surface. This is condition for long-lasting plastic elasticity. Eartips bodies in case of washing should be unscrewed. The membrane can be removed after carefully pushing away the plastic ring and washed in soapy water. A replacement membrane is supplied. All components can be disinfected with ethyl or isopropyl alcohol. Do not use thermal sterilization.

#### Warranty

Manufacturer warranty: for cardiology stethoscopes for a period of 5 years, for the rest of stethoscopes for 3 years from the date of sale. We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. This guarantee excludes any defects which occur because of normal wear and tear, incorrect use, the use of non-original spare parts, or non-authorized modification or repairs. For repairs, please contact your supplier, who will return the goods to us.

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